

Safety Information

Table of contents

1	Safety Information	1
1.1	Important Safety Instructions	1
1.2	General Warnings	1
1.2.1	Electrical Safety	1
1.2.2	Patient Safety	1
1.2.3	Optical Safety	2
1.3	Infection Control	2
1.3.1	Required Procedures	2
1.3.2	Approved Disinfectants	2
1.4	Environmental Conditions	2
1.4.1	Operating Environment	2
1.4.2	Storage Environment	2
1.5	Disposal	2
1.5.1	End of Life	2
1.5.2	Battery Disposal	2
1.6	Reporting	2

1 Safety Information

1.1 Important Safety Instructions

! Critical Safety Information

Read and understand all safety information before using the OralEye system. Failure to follow these instructions may result in injury or device damage.

1.2 General Warnings

1.2.1 Electrical Safety

- Use only the provided power adapter
- Do not expose the device to liquids
- Disconnect power before cleaning
- Do not use if power cord is damaged

1.2.2 Patient Safety

- Ensure patient comfort during imaging procedures
- Avoid prolonged LED exposure to eyes
- Use appropriate infection control procedures
- Follow standard precautions for patient contact

1.2.3 Optical Safety

- LED light source is safe for intended use
- Do not stare directly into LED array
- Automatic exposure control prevents overexposure

1.3 Infection Control

1.3.1 Required Procedures

1. **Barrier Protection:** Use disposable barriers for intraoral components
2. **Surface Disinfection:** Clean all surfaces with approved disinfectant after each patient
3. **Hand Hygiene:** Follow standard hand hygiene protocols

1.3.2 Approved Disinfectants

- Isopropyl alcohol (70%)
- EPA-registered hospital-grade disinfectants
- Follow manufacturer's contact time recommendations



Warning

Do not autoclave or immerse device in liquid. Do not use abrasive cleaners.

1.4 Environmental Conditions

1.4.1 Operating Environment

- **Temperature:** 15°C to 30°C (59°F to 86°F)
- **Humidity:** 30% to 75% relative humidity
- **Pressure:** 700 hPa to 1060 hPa

1.4.2 Storage Environment

- **Temperature:** -10°C to 50°C (14°F to 122°F)
- **Humidity:** 10% to 90% relative humidity (non-condensing)

1.5 Disposal

1.5.1 End of Life

- Do not dispose of device in regular waste
- Follow local regulations for electronic waste disposal
- Contact manufacturer for recycling options

1.5.2 Battery Disposal

- Lithium-ion batteries must be recycled according to local regulations
- Contact local waste management for battery recycling

1.6 Reporting

Report adverse events or device problems to:

- **Manufacturer:** support@oraleye.ai
- **FDA:** via MedWatch (1-800-FDA-1088)

This safety information is current as of the document effective date.